

VRN'S NEW "SONIC" ERA SCALER

POWER

6

0

5



Built-in High-Performance Chip

Equipped with an STMicroelectronics chip, featuring advanced algorithms that automatically detect resonance frequency, preventing instantaneous impact and ensuring a comfortable scaling experience.



High-Quality Lighted Handpiece

New handpiece assembled technology with a detachable design, able to use high-temperature and high-pressure sterilization. Soft light source provides clearer visibility during treatment.



Smart Touch System

Responsive waterproof touch panel with one-touch switching between **G** (scaling), **P** (periodontal), and **E** (endodontic) modes. Offers **12 power** levels for precise treatment customization.



Automatically adjusts power based on calculus hardness, increasing power as needed for enhanced cleaning efficiency.

SUPRAGINGIVAL SCALING



Before







Supragingival Scaling
G1
Supragingival scaling, used to removing in all directions.
G2
Flat edge is used to remove heavy su
G4
Used to remove all supragingival cale and supragingival interdental areas.
Periodontal Treatment
P1 Used for subgingival calculus remova



Manufacturer Name: Guilin Veirun Medical Technology Co., Ltd. Address: No.D-07 Information Industry District, High-Tech Zone, 541004 Guilin, Guangxi, PEOPLE'S REPUBLIC OF CHINA **Tel:** +86(773)5880619 **Tel:** +86(773)2260566 Fax: +86(773)2260519 Web: www.veirun.com E-mail: info@veirun.com





Scan to visit VRN official website



VEIRUN.COM