

Ultrasonic Periodontal Therapy System

FIVE-IN-ONE EXCELLENT INTEGRATION


- Ultrasonic scaling
- Periodontal therapy
- Root canal irrigation
- Supragingival airpolishing
- Implant maintenance.



DUAL-DISC MULTI-POINT INTELLIGENT TOUCH PANEL

Fully functional touch screen with 12 power and water levels, offering a wide adjustable range. G Liquid NBX Min BOMEL


INTELLIGENT MULTI-CORE CHIP



ELLIPTICAL AMPLITUDE, COMFORTABLE SCALING

000

Elliptical vibration trajectory enables 0° angle operation of the scaling tip on the tooth surface, providing polishing effect, dispersing the impact force between the scaling tip and the tooth, achieving comfortable and 360° dead angle-free scaling.




AUTOMATIC IDENTIFICATION OF CLEANING MODE

000

The system can automatically identify the scaling or airpolishing mode based on the selected working handpiece.





HIGH-PRECISION MECHANICAL PRESSURE ADJUSTMENT

0 0 0


Accurately controls the air pressure range, delivering stable output, efficient operation, and a comfortable experience.



CONSTANT PRESSURE POWDER ADJUSTMENT

0 0 0

Easy control, comfortable and efficient. Dentists can adjust the output of powder as needed under fixed pressure, reducing discomfort caused by air pressure changes.




Ordinary airpolishing machines change the powder output by adjusting the air pressure, which causes uncomfortable such as tooth soreness during cleaning.

With DQ-40, the powder output is adjustable under fixed air pressure.


DQ-40 and other sand blasting machine sand output comparison

Adjust the amount of sand at will

Knob lid for easy handling








HIGH-PRESSURE REVERSE BLOWING


0 0 0

Patented technology, no need to use sand cleaning needle for handpiece cleaning, quickly removes sand powder accumulation.



EFFICIENT REMOVAL OF DENTAL PLAQUE

0.7mm aperture supragingival airpolishing head, precision-made with imported German steel, reliable and durable, more concentrated airpolishing, higher efficiency and more comfortable teeth cleaning.



AUTOMATED DUAL PIPELINE CLEANING

Newly upgraded airpolishing pipeline one-click cleaning, efficiently removes residual prophylaxis



TITANIUM ALLOY SCALING TIPS

0 0 0


Equipped with titanium alloy periodontal tips, offering low hardness to ensure no harm to the teeth.



MULTIFUNCTIONAL WIRELESS FOOT PEDAL

0 0

Wireless era, unlimited possibilities, eliminating cable tangles.



IPX4



Manufacturer Name: Guilin Veirun Medical Technology Co., Ltd. Address: No.D-07 Information Industry District, High-Tech Zone,

541004 Guilin, Guangxi, PEOPLE'S REPUBLIC OF CHINA Tel: +86(773)5880619

Tel: +86(773)2260566
Fax: +86(773)2260519
Web: www.veirun.com
E-mail: info@veirun.com

Scan to visit VRN official website



VEIRUN.COM